



**Bledsoe  
Brace  
Systems™**

## Bledsoe Hi-Top Boot



**Superior support without  
casting complications**

- **Reformable:** Aluminum boot shell and upright system that unlike plastic boots, can be reformed to fit virtually any leg or ankle shape.
- **Low Profile:** Low profile rocker bottom helps eliminate pelvic tilt and alleviate low back pain.
- **Patient comfort:** Ultra breathable foam wrap for maximum comfort.
- **Fast, easy fitting:** Saves time over complicated casting procedures.
- **Unique rocker design:** Allows patients to closely simulate a near normal gait without abnormal loads on any joints.
- **Hinged Boot:** Hinged Hi-Top boot is available that provides R.O.M. from 20° dorsi-flexion to 40° plantar-flexion in 10° increments.
- **Optional Air Ankle / Heel Pad:** Allows for ability to adjust ankle pressure for increased comfort and provides pneumatic ankle support.



# Bledsoe Hi-Top Boot

## Restrictions on use

This device is offered for sale by or on the order of a physician or other qualified medical professional. This device is not intended for reuse on a second patient or for extended use beyond 180 days. This device is generally not intended for patients in excess of 250 lbs. (114 kg.) This device is not intended to prevent injuries or reduce or eliminate the risk of injury.

## Warning

Do not leave this device in the trunk of a car on a hot day or in any place where the temperatures may exceed 150°F (65°C) for any extended period of time. If any additional pain or symptoms occur while using this device, discontinue use and seek medical attention.



**Bledsoe  
Brace  
Systems™**

by Medical Technology, Inc.

2601 Pinewood Drive  
Grand Prairie, TX 75051 USA

1-888-BLEDSOE  
(1-888-253-3763)

1-800-527-3666

972-647-0884

Local Fax:  
972-660-5495

International Fax:  
972-606-0649

[www.bledsoebrace.com](http://www.bledsoebrace.com)

## Indications

Indicated for post-injury or post-operative application to grade 1, 2, or 3 ankle sprains, and stable or internally fixed fractures of the foot or ankle. Also indicated for stable or internally fixed fractures of the middle to distal tibia or fibula when the optional tibial fracture bracing plate is used.

## Contraindications

Contraindicated for unstable fractures and proximal tibial or fibular fractures.

## Bledsoe Hi-Top Boot vs. other methods

Features	Bledsoe Boots	All other low profile walkers	All other profile walkers
Unique low profile rocker design allowing for a near normal walking gait and alleviating lower back pain.	Yes	No (rocker too flat)	No (rocker curved, but too high)
Aluminum shell and uprights permit forming to fit most any size ankle or leg.	Yes	No	No
Ultra breathable foam wraps optimize patient comfort and skin protection.	Yes	No	No
Open back design allows adjustment of foot and leg position over rocker to provide best walking gait.	Yes	No	No

## Ordering Information Hi-Top Boot

Shoe size (preferred) or leg measurement taken from knee joint line to medial malleolus.

Shoe Size			Leg Measurement	Boot	Fixed Hi-Top		EZ Set Hi-Top	
Men	Women	Euro	Length	Size	W/Cuff	W/Bootie	W/Cuff	W/Bootie
1- 2.5	1-3.5	29-33	9"-11" (23 - 28 cm)	XS	AL?36001	AL?36021	AL?72001	AL?72021
2.5-5.5	3.5-6	33-37	11"-13" (28 - 33 cm)	S	AL?36003	AL?36023	AL?72003	AL?72023
5.5-9	6-9.5	37-41	13"-15" (33 - 38 cm)	M	AL?36005	AL?36025	AL?72005	AL?72025
9-12.5	9.5-13	41-46	15"-17" (38 - 43 cm)	L	AL?36007	AL?36027	AL?72007	AL?72027
12.5-17	13 +	46-49 +	17"-19" (43 - 48 cm)	XL	AL?36009	AL?36029	AL?72009	AL?72029

**Note: When ordering replace the ? with one of the following: 0 = ankle / heel pad 2 = air ankle / heel pad**

	Hinge & Thigh Attachment			Hinge Choice:		
	If Leg Length Is	Use Size	Simple	Simple w/ QuikLok	Multi-Centric	
This product consists of thigh portion and hinges of the Bledsoe Brace but without calf components. It is designed to attach to the Hi-Top Boot to create a long leg device.	17" - 19" (43 - 48 cm)	18" (46 cm)	AK025018	AK0525018	AK015018	
	19" - 21" (48 - 53 cm)	20" (51 cm)	AK025020	AK0525020	AK015020	
	21" - 23" (53 - 59 cm)	22" (56 cm)	AK025022	AK0525022	AK015022	
	23" - 25" (59 - 64 cm)	24" (61 cm)	AK025024	AK0525024	AK015024	
	25" - 27" (64 - 69 cm)	26" (66 cm)	AK025026	AK0525026	AK015026	
	27" - 29" (69 - 74 cm)	28" (71 cm)	AK025028	AK0524028	AK015028	
	29" - 31" (74 - 79 cm)	30" (76 cm)	AK025030	AK0525030	AK015030	
	31" - 33" (79 - 84 cm)	32" (81 cm)	AK025032	AK0525032	AK015032	
33" - 35" (84 - 89 cm)	34" (86 cm)	AK025034	AK0525034	AK015034		

**Note: Leg measurement for Hinge & Thigh Attachment is taken from 1" (2.5cm) below groin on medial side of leg to peak of medial malleolus.**

**Accessories:** Rain, Cold Weather, or Night Cover - Tibial Plates - Replacement cuff and booties, Strap kits  
**Call for more information.**

## Bledsoe limited warranty

This device is warranted to be free from defects in material and workmanship for a period of 90 days. These warranties apply to devices that have not been modified or subjected to misuse, abuse or neglect.



For more detailed information on this brace including fitting and adjustment instructions, visit our Web site: [www.bledsoebrace.com](http://www.bledsoebrace.com)

**For product information or questions pertaining to sales or service, please contact the sales representative in your area or Bledsoe Brace Systems.**

