



**Bledsoe
Brace
Systems™**

Bledsoe Simple Pelvic Brace



Sturdy,
formable
aluminum
uprights

Soft
breathable
foam
wraps

- **Rotational control:** Controls rotation, abduction and limits or locks range of motion about the hip.
- **Shape to fit:** Sturdy, formable aluminum uprights and rigid plastic shells for positive control.
- **Patient comfort:** Soft breathable foam wraps for patient comfort. Non-stretch pile straps that encompass the entire brace, interlocking all components.
- **The Bledsoe Pelvic Brace:** Limits posterior flexion to 30° and anterior flexion to 105° or locks motion from 30° posterior flexion to 60° anterior flexion in 15° increments.
- **Short Application:** This off-the-shelf brace fits like a custom and can be applied in less than 10 minutes without tools.

**Solves the problem of hip
rotational control**

Bledsoe Simple Pelvic Brace

Restrictions on use

This device is offered for sale by or on the order of a physician or other qualified medical professional. This device is not intended for reuse on a second patient. This device is not intended to prevent injuries or reduce or eliminate the risk of injury.

Warning

Do not leave this device in the trunk of a car on a hot day or in any place where the temperatures may exceed 150°F (65°C) for any extended period of time. If any additional pain or symptoms occur while using this device, discontinue use and seek medical attention.



**Bledsoe
Brace
Systems™**

by Medical Technology, Inc.

2601 Pinewood Drive
Grand Prairie, TX 75051 USA

1-888-BLEDSOE
(1-888-253-3763)

1-800-527-3666

972-647-0884

Local Fax:
972-660-5495

International Fax:
972-606-0649

www.bledsoebrace.com

Indications

Indicated to limit motion and reduce incidence of dislocation of the hip and to provide limited motion of hip abduction or adduction after injury or operative procedures to the hip joint.

Contraindications

Contraindicated for unstable fractures of the pelvis or femur.

Ordering Information

Simple Pelvic Brace

Measure the leg length and use the conversion chart shown below to determine Pelvic Brace size. If the iliac crest or thigh circumference is larger than maximum, you must order a custom brace. **Do not order the next larger size!** Example: If the patient has a 22" (56 cm) leg length, but their iliac crest circumference is 43" (109 cm), they will need a custom Pelvic Brace.

<p>To measure leg length: Measure from crotch to ankle as shown.</p>	Leg Length	Iliac Crest Circumference	Thigh Circumference	Size	Left	Right
	(56 - 66 cm) 22" - 26"	(96.5 cm Max.) 38" Max.	(66 cm Max.) 26" Max.	S	AT013101	AT013201
	(66 - 76 cm) 26" - 30"	(114 cm Max.) 45" Max.	(76 cm Max.) 30" Max.	M	AT013103	AT013203
	(76 - 86 cm) 30" - 34"	(140 cm Max.) 55" Max.	(86 cm Max.) 34" Max.	L	AT013105	AT013205

Accessories

• Abdominal cuff and strap kits

• Thigh cuff and strap kits

Bledsoe limited warranty

This device is warranted to be free from defects in material and workmanship for a period of 90 days for wraps, straps and normal wear components, 12 months on all other parts. These warranties apply to devices that have not been modified or subjected to misuse, abuse or neglect.



For more detailed information on this brace, including fitting and adjustment instructions, visit our Web site:

www.bledsoebrace.com

For product information or questions pertaining to sales or service, please contact the sales representative in your area or Bledsoe Brace Systems.



CP010117--- Rev.G 6/04